Material included in the package:

Product code: 4 SCP-20

1 x 10 ml GM501 SpermStore

Material not included in the package:

Sperm freezing straws

Sterile pipettes

Freezing tank with liquid nitrogen

LAF bench (ISO 5 environment)

Sperm wash medium

Product specifications and quality control:

- All raw materials are of highest available purity (European Pharmacopoeia and/or USP standard) if applicable.
- A certificate of analysis is issued for each batch upon request from our supplier with respective lot number.
- The MSDS for GM501 SpermStore is available upon request and can also be downloaded from our website.

GM501 SpermStore is manufactured and tested according to the following specifications:

pH	7.20-7.90 (Release criteria: 7.20-7.60)
Sterility	sterile - SAL10 ⁻³ (Sterility Assurance Level)
Endotoxins (EU/ml)	< 0.25

Sperm Survival test

>=0 survival after 4 hours exposure of density selected spermatozoa to the test medium in %

Mouse Embryo Assay (MEA)

not MEA tested

EN: Intended use/Intended users:

- GM501 SpermStore is intended for freezing human sperm including epididymal and testicular sperm.
- The intended users are IVF professionals (lab technicians, embryologists or medical doctors).

Composition:

GM501 SpermStore is a ready-to-use HEPES buffered cryopreservation medium which also contains physiologic salts, glycine, dextrose, monohydrate, lactate, glycerol, sucrose and human serum albumin (4.00 g/liter) to protect sperm from damage due to the freezing procedure.

GM501 SpermStore and Sperm Preparation:

GM501 SpermStore can be used in combination with GM501 Gradient.

Instructions for use:

- Ensure all media are well mixed before use.

Before freezing:

- In case of very low sperm concentrations it is advisable to concentrate the sperm before freezing. In case of very high sperm concentration, sperm can be diluted using a suitable sperm washing medium.
- GM501 SpermStore needs to be warmed to room temperature before use to avoid cold-shock.

Freezing:

- Allow the semen to liquefy at room temperature for 30 minutes.
- Mix 1 ml of sperm with 0.7 ml of GM501 SpermStore. Add the medium in drops while gently swirling.
- Leave the mixture for 10 minutes at room temperature for equilibration.

Aspirate the sample/medium mixture into the freezing straws, leaving approximately 1.5 cm of air at the end of the straws.

Seal the straws.

Dry off individually with a lint-free wipe.

Shake the straw to the centre of the straw.

Place the straws horizontally on a styrofoam board in a liquid nitrogen bath to allow for freezing in vapour phase. Leave for (at least) 15 minutes.

Transfer straws quickly into liquid nitrogen and store at -196°C.

Thawing:

- Remove as many straws as required from the liquid nitrogen.

Place the samples in tap water for 5 minutes (room temperature or 37°C).

Cut off one end of straw, place the open end inside a container (e.g. a test tube) and tap the straw lightly, against the side of the container to allow complete evacuation of the mixture.

Dilute the thawed sperm in a suitable insemination medium (at least 3 ml per 0.5 ml semen) and mix thoroughly.

Centrifuge for 15 minutes at 300-350 x g.

Resuspend pellet in a suitable insemination medium.

After thawing:

- If necessary, use sperm preparation techniques after thawing the semen to eliminate dead sperm cells and debris. Dilute the concentrated sperm in a suitable insemination medium.

Precautions and warnings:

- Standard measures to prevent infections resulting from the use of medical products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medical products are prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot totally be excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes.

Therefore, handle all specimens as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling specimens.

Product does not contain any antibiotics. Always work under strict hygienic conditions (e.g. LAF-bench ISO Class 5) to avoid potential contamination.

Only for the intended use.

The user of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

Pre-use checks:

- Do not use the product if it becomes discoloured, cloudy, or shows any evidence of microbial contamination.
- Do not use the product if seal of the container is opened or defective when the product is delivered.

Storage instructions and stability:

- The shelf life is 18 months from time of manufacture.
- Store at 2-8°C.
- Do not freeze before use.
- Keep away from (sun) light.

The product can be used safely up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8°C.

Stable after transport (max. 5 days) at elevated temperature (≤ 37°C).

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