

# Produktbeschreibung

## Use less toxin & Exact Dosing

Low dead space save up to 20%. Predictable and reproducible, no over- and under-doses.

## Focus

Tactile and audible feedback. Focus only on patient and case.

## Calculation ml-unit

Is included in the syringe and setting.



## Locking handle for

- Free moving-mode: plunger move free (pull up);
- Click-mode: operational, one click per dose.

## Ergonomic handle

- Handle;
- Thumb-rest.

## Display

- Remaining units;
- Stand upright.

